

<p>Site de Fribourg: Microbiologie Médicale et Moléculaire, Université de Fribourg Chemin du Musée 18 1700 Fribourg Tel: +41 26 300 9581 e-mail : nara@unifr.ch</p>	 <h1 style="margin: 0;">NARA</h1> <p style="margin: 0;">Nationales Referenzlaboratorium zur Früherkennung neuer Antibiotikaresistenzen und Resistenzmechanismen</p>	<p>Site de Lausanne: Laboratoire d'épidémiologie CHUV Rue du Bugnon 46 1011 Lausanne Tel: +41 21 314 02 60 email: dominique.blanc@chuv.ch</p>  <p style="text-align: right; font-size: small;">ISO 17025:2005</p>
--	---	---

10/05/2019

Polymyxin resistance in Gram negatives: identification and screening

Prof. Patrice Nordmann, Dr. Dominique Blanc, Dr. Laurent Poirel

Multidrug resistant (MDR) bacteria are of global concern, notably with the spread of carbapenemase producers. This is one of the reason why polymyxins such as colistin, that are old antibiotics, regain popularity as a last resort treatment to face the worldwide emergence of MDR Gram-negative bacteria. Colistin is a polycationic and bactericidal drug that targets the lipid A moiety of the lipopolysaccharride (LPS) moving its cationic charges and subsequently leading to cell wall lysis and bacterial death. The increasing use of colistin in animals since years but also nowadays in humans has led to emerging resistance to that antibiotic. Until recently, all mechanisms described involved chromosomal genes encoding proteins mainly related to regulatory systems of the LPS biosynthesis, such as the PmrAB and PhoPQ two-components systems in Enterobacteriales. In November 2015, the first plasmid-mediated colistin resistance gene, namely *mcr-1* that encodes a phosphoethanolamine transferase, was identified in enteobacterial isolates. Then *mcr*-like genes, namely *mcr-2* to *mcr-8*, have been identified. The *mcr* genes have been reported worldwide mostly in animal where the selection pressure by using polymyxin is high, and among environmental isolates. Although they are reported in clinical strains from humans in Switzerland, their spread among human pathogens in this country remains rare according to the results of recent surveys.

Resistance to colistin in human *E. coli* is mostly related to non-MCR and chromosome-encoded mechanisms. We are currently observing a worldwide spread of colistin-resistant strains (with a chromosome-encoded mechanism of resistance) among carbapenemase producers (KPC, OXA-48, NDM). Such strains have been identified in Switzerland resulting from an importation process (mostly from Italy) but also from direct selection in Switzerland following the treatment of patients with colistin.

A particular concern is that the majority of the MCR producers may be difficult to detect since they confer low MIC values (around 4 mg/L) while chromosomally-encoded resistance mechanisms confer much higher levels of resistance (usually > 64 mg/L).

So far resistance to polymyxins in *Pseudomonas aeruginosa* and *Acinetobacter baumannii*, is only chromosome-encoded with the exception of few plasmid-encoded MCR-4 in that latter species. They are mostly related to changes in LPS structure.

A-Infection

Testing shall be performed from isolated colonies

1- Identification

- MALDI-TOF
- or at the species level: Vitek (bioMérieux) or Phoenix (PMIC/ID Panels) for example

2- Susceptibility testing

2.1 MIC testing remains the gold standard for determining resistance to polymyxins. It must always be determined in liquid medium in cation-adjusted with a final composition of 20-25 mg/L of calcium and 10-12.5 mg/L of magnesium. Colistin sulphate must be used in untreated polystyrene plates without addition of any surfactant. One shall pay attention to the purity of the colistin (a natural product) that varies from one furnisher to another. We are using colistin tablets (ref. TAB/CO 0,8, Mast Group) as primary source of antibiotic. The E-test technique is not recommended for detecting colistin resistance due to the poor diffusion of colistin in solid agar.

The MIC breakpoints of colistin for Enterobacteriales, *Pseudomonas aeruginosa* and *Acinetobacter baumannii* are indicated below (EUCAST 2019 guidelines v.9.0)/

:

	MIC breakpoints (mg/L)	
	S ≤	R >
Colistin	2	2


Results of susceptibility testing are obtained in 18 h. Since the MCR proteins are zinc-dependant enzymes, several studies reported the possibility to inhibit MCR activity and therefore to decrease the MICs of colistin by adding dipicolinic acid and EDTA. This inhibition may work for detecting MCR producers but only in *E. coli* and only with dipicolinic acid.

2.2. MIC determination is proposed by several companies. Those commercial devices are UMIC Colistine (Biocentric, France), MIC Strip Colistin (MERLIN, Germany), SensiTest Colistin (Liofilchem, Italy), Sensititre system (Thermo Fisher Scientific USA) and some automated systems. The UMIC Colistine consists of unitary tests containing dehydrated concentrations of colistin. It has been reported as reliable. The MIC strip colistin is based on the same principle. The Sensitest Colistin based also on the same principle consists of compact panel of 4 tests. The Sensititre system includes different antibiotics on 96-well trays with a customizable plate layout. The automated systems were developed to shorten results timeframes by increasing sensitivity, and also to avoid manipulation bias with incubation and real-time reading. However, they cannot give a real MIC value. The manufacturers are MicroScan Walkaway (Beckman Coulter, USA), Vitek (bioMérieux, France) and BD Phoenix (Becton Dickinson, France). Their performance are variable in particular for detecting resistance to polymyxins in non fermenters. With the exception of the Sensitest and the Sensititre methods, very major errors may occur with reports of false sensitivity for MCR producers, those latter often exhibiting low-level colistin resistance.

2.3. Rapid susceptibility testing; qualitative detection of colistin resistance

Rapid susceptibility tests for detection of colistin resistance have been developed for Enterobacteriales, *P. aeruginosa* and *A. baumannii* to provide faster results, obtained in 2-4 h. The Rapid Polymyxin NP test is based on a change of pH according to glucose metabolism related to bacterial growth. Detection of colistin resistance is made according to growth in presence of a given concentration of colistin. The test has been developed from bacterial cultures and can also be used directly using positive blood cultures. Compared to the broth microdilution susceptibility method, there is excellent agreement for detecting MCR-1 and MCR-2 producers. The Rapid Polymyxin NP test has also a good sensitivity to detect naturally-occurring polymyxin resistance mechanisms such as those observed in *Hafnia alvei*. However it fails to detect several colistin-resistant *Enterobacter* sp. isolates probably due to the frequent colistin heteroresistance pattern observed in this species. This test has been commercialized under the trade name of Rapid Polymyxin NP test (ELITECH Microbiology).

A rapid test was developed also for identification of polymyxin resistance in non fermenters (Ras Polymyxin *Acinetobacter/Pseudomonas* NP). This test detects viable cells after growth in a medium containing a defined concentration of colistin. The principle of this test is based on the visual detection of the reduction of the resazurin reagent, a viability colorant, as observed by its color change (blue to purple or pink). Its evaluation was performed by using colistin-resistant and colistin-susceptible *A. baumannii* and *P. aeruginosa* isolates. Sensitivity and specificity were found to be 100% and 95%, respectively, by comparison with the standard broth microdilution method. The Rapid Resa Polymyxin *Acinetobacter/Pseudomonas* NP test is inexpensive, easy to perform, highly sensitive and specific, and can be completed in 4 h. A derivative of the Rapid Polymyxin NP test has been developed (ELITECH Microbiology) that also shows good performances for detecting colistin resistance in those species.






Figure 1. Rapid Polymyxin NP for detecting colistin resistance in Enterobacteriales. Note the lack of growth of colistin-susceptible strain in presence of colistin.

A MALDI-TOF technique (MALDIxin test) has been developed to detect chromosome and plasmid-mediated resistance in Enterobacteriaceae. It detects this resistance trait in 15 min. However possessing a MALDI-TOF machine is mandatory and the technique has to be settled in the lab.

- 3- Immunological detection of MCR producers. Lateral flow techniques have been commercialized for detecting MCR-1 producers such as that from NG Biotech (France). The test is rapid (15 min), sensitive and specific. It does not need any additional equipment for the reading.



Figure 2. Lateral-flow techniques for detecting MCR-1 producer

- 4- Molecular identification of polymyxin resistance. These methods are complementary to the phenotypic techniques. The main mutations for enterobacteriales are located on genes coding the two-component systems PmrA/PmrB and PhoP/PhoQ. Mutations in the negative feedback regulator of PhoPQ is the main mechanism of resistance to polymyxins in *K. pneumoniae*. Gene amplification and sequencing takes usually at least 3 days. Primers have been developed to identify the *mcr* genes and multiplex PCR genes for identification of those genes.

Real-time PCR techniques have been also developed for detection of *mcr-1/mcr-2* genes with excellent sensibility and specificity. Semi-automated PCR-based detection of *mcr* genes have been developed such as the CT103XL® microarray and the Eazyplex® (Amplex Biosystems, Germany) with good sensibility and specificity but remained not cost-effective.

Whole genome sequencing (WGS) may also be used, in particular for detecting unknown mechanisms of resistance to polymyxins. WGS so far can only be performed by expert laboratories and remains neither cost effective nor rapid.

We recommend therefore the following strategy for detection of polymyxin resistance;

- 1- Bacterial identification
- 2- Rapid Polymyxin NP test (and in the special case of detection of MCR- only, lateral flow technique)
- 3- Confirmation of resistance by performing susceptibility testing
- 4- PCR and sequencing of *mcr* genes or defined chromosomal genes to precisely determine the molecular bases of resistance (as determined by the NARA).

B - Screening (colonization)

- 1- Which patients ?

No recommendation consensus has been established for screening patients carrying colistin-resistant Gram-negative isolates. However, we recommend to screen contact patients of colonized or infected patient with a colistin- and carbapenem-resistant enterobacterial strain, or colistin-resistant and multidrug-resistant *P. aeruginosa* or *A. baumannii*.

- 2- Which samples ?

- Stools
 - Rectal swabs of good quality
- Repeated screening of samples may be recommended if patients are treated with antibiotics (false-negatives results) and for high-risk patients.*

- 3- Which screening medium ?

- Several screening media have been developed for screening. All of them contain at least a selective concentration of colistin. The first screening medium that had been developed is the SuperPolymyxin medium subsequently commercialized by ELITECH Microbiology. It has the advantage of facilitating the visualization of *E. coli* strains (metallic green reflect) because its composition (EMB agar). It detects fermenters and

non-fermenters and strains displaying low-level resistance colistin. The CHROM agar COL-ASPE medium present the advantage to be chromogenic with the capacity to differentiate colistin-resistant non fermentative Gram-negative strains as well as Enterobacteriales. Recently, bioMérieux has also commercialized a screening medium. CHROMID R colistin R which detects colistin-resistant enterobacteriales.

- We showed very recently that recovery of colistin-resistant enterobacteriales may be enhanced by a pre-culture step (18 h) of stools in a broth containing 1 mg/L of colistin (unpublished data).

4- Recommendations

- Non outbreak situation (Fig 3); direct plating of rectal swabs or stools on chromogenic selective media.
- Outbreak situation (Fig. 4)
 - Direct plating on chromogenic selective media and/or use of molecular techniques if available
 - Concomitant cultures of the samples in an enrichment broth containing 1 mg/L of colistin.




Figure 3. Detection of carriage of colistin-resistant Gram negatives (in the absence of an outbreak context)




Figure 4. Detection of carriage of colistin-resistant Gram negatives (in the presence of an outbreak context)

Recommended reading

1. Behera B, Mathur P, Das A, Kapil A, Gupta B, Bhoi S, Farooque K, Sharma V, Misra MC. 2010. Evaluation of susceptibility testing methods for polymyxin. *Int J Infect Dis* 14:e596 – e601.
2. Bernasconi O, Principe L, Tinguely R, Karczmarek A, Perreten V, Endimiani A. 2017. Evaluation of a new commercial microarray platform for the simultaneous detection of β -lactamase and *mcr-1* and *mcr-2* genes in *Enterobacteriaceae*. *J Clin Microbiol.* 55:3138-3140
3. Biswas S, Brunel JM, Dubus JC, Reynaud-Gaubert M, Rolain JM. 2012. Colistin: an update on the antibiotic of the 21st century. *Expert Rev Anti Infect Ther* 10:917–934
4. Bontron S, Poirel L, Nordmann P. 2016. Real-time PCR for detection of plasmid-mediated polymyxin resistance (*mcr-1*) from cultured bacteria and stools. *J. Antimicrob. Chemother.* 71:2318-2320.
5. Bourrel A.S., Poirel L., Royer G., Darty M., Vuillemin X., Kieffer N., Clermont O., Denamur E., Nordmann P., Decousser J.W. (2019). Colistin resistance in Parisian inpatient faecal *Escherichia coli* as the result of two distinct pathways. *J Antimicrob Chemother.* doi 10.1093/jac/dkz090.
6. Büdel T, Clement M, Bernasconi OJ, Principe L, Perreten V, Luzzarro F, Endimiani A. 2018. Evaluation of EDTA- and DPA-based microdilution phenotypic tests for the detection of MCR-mediated resistance to colistin in *Enterobacteriaceae*. *Microb Drug Resist* doi:10.1089/mdr.2018.0275.
7. Cannatelli A, Di Pilato V, Giani T, Arena F, Ambretti S, Gaibani P, D'Andrea MM, Rossolini GM. 2014. In vivo evolution to colistin resistance by PmrB sensor kinase mutation in KPC-producing *Klebsiella pneumoniae* is associated with low-dosage colistin treatment. *Antimicrob Agents Chemother* 58:4399 – 4403.

8. Cannatelli A, Giani T, D'Andrea MM, Di Pilato V, Arena F, Conte V, Tryfinopoulou K, Vatopoulos A, Rossolini G. 2014. MgrB inactivation is a common mechanism of colistin resistance in KPC-producing *Klebsiella pneumoniae* of clinical origin. *Antimicrob Agents Chemother* 58:5696–5703.
9. Carretto E, Brovarone F, Russello G, Nardini P, El-Bouseary M, Ablokaish AF, Tyrell JM. 2018. Clinical validation of the SensiTest colistin, a broth microdilution based method to evaluate colistin MICs. *J Clin Microbiol* 26:doi:10.1128/JCM.01523-17.
10. Cheng YH, Lin TL, Pan YJ, Wang YP, Lin YT, Wang JT. 2015. Colistin resistance mechanisms in *Klebsiella pneumoniae* strains from Taiwan. *Antimicrob Agents Chemother* 59:2909–2913.
11. Dona V, Bernasconi O, Kaisraian S, Tinguely R, Endimiani R. 2017. A SYBR® green based real-time PCR methods for improved detection of mcr-1 mediated colistin resistance in human stool samples. *J Glob Antimicrob Resist* 9:57-60.
12. Falagas ME, Rafailidis PI, Matthaïou DK. 2010. Resistance to polymyxins: mechanisms, frequency and treatment options. *Drug Resist Updat* 13:132–138.
13. Giani T, Arena F, Vaggelli G, Conte V, Chiarelli A, Henrici De Angelis L, Fornaini R, Grazzini M, Niccolini F, Pecile P, Rossolini GM. 2015. Large nosocomial outbreak of colistin-resistant, carbapenemase-producing *Klebsiella pneumoniae* traced to clonal expansion of an *mcrB* deletion mutant. *J Clin Microbiol* 53:3341–3344.
14. Imirzalioglu C, Falgenhauer L, Schmiedel J, Waezsada SE, Gwodzinski K, Roschanski N, Chakraborty T. 2017. Evaluation of a loop-mediated isothermal amplification-based assay for the rapid detection of plasmid-encoded colistin resistance gene *mcr-1* in enterobacteriaceae isolates. *Antimicrob Agents Chemother* 61;poee/02326-16.
15. Jayol A, Dubois V, Poirel L, Nordmann P. 2016. Rapid detection of polymyxin-resistant *Enterobacteriaceae* from blood cultures. *J Clin Microbiol* 54:2273–2277.
16. Jayol A, Kieffer N, Poirel L, Guérin F, Güneser D, Cattoir V, Nordmann P. 2018. Evaluation of the Rapid Polymyxin NP test and its industrial version for the detection of polymyxin-resistant *Enterobacteriaceae*. *Diagn Microbiol Infect Dis*,doi:10.1016.
17. Jayol A, Nordmann P, Lehours P, Poirel L, Dubois V. 2017. Comparison of methods for detection of plasmid-mediated and chromosomally encoded colistin resistance in *Enterobacteriaceae*. *Clin Microb Infect* 24/30291-30294.
18. Jayol A, Nordmann P, André C, Poirel L, Dubois V (2018). Evaluation of three broth microdilution systems to determine colistin susceptibility of Gram-negative bacilli. *J Antimicrob Chemother* 73:1272-1278.
19. Jayol A, Poirel L, Brink A, Villegas MV, Yilmaz M, Nordmann P. 2014. Resistance to colistin associated with a single amino acid change in protein PmrB among *Klebsiella pneumoniae* isolates of worldwide origin. *Antimicrob Agents Chemother* 58:4762–4766.
20. Jayol A, Dortet L, Poirel L, Nordmann P. 2016. National survey of colistin resistance among carbapenemase-producing *Enterobacteriaceae* and outbreak caused by colistin-resistant OXA-48-producing *Klebsiella pneumoniae*, France. 2014. *Euro Surveill* doi:10.2807/1560-7917.
21. Jayol A, Poirel L, Villegas MV, Nordmann P. 2015. Modulation of *mcrB* gene expression as a source of colistin resistance in *Klebsiella oxytoca*. *Int J Antimicrob Agents* 46:108–110.
22. Humphries RM. 2015. Susceptibility testing of the polymyxins: where are we now? *Pharmacotherapy* 35:22–27.
23. Lee SY, Shin JH, Lee K, Joo MY, Park KH, Shin MG, Suh SP, Ryang DW, Kim SH. 2013. Comparison of the Vitek 2, MicroScan, and Etest methods with the agar

- dilution method in assessing colistin susceptibility of bloodstream isolates of *Acinetobacter* species from a Korean university hospital. *J Clin Microbiol* 51:1924 – 1926.
24. Landman D, Georgescu C, Martin DA, Quale J. 2008. Polymyxins revisited. *Clin Microbiol Rev* 21:449 – 465.
 25. Lescat M, Poirel L, Jayol A, Nordmann P. 2018. Performances for the rapid polymyxin *Acinetobacter* and *Pseudomonas* tests for colistin susceptibility testing. *Microb Drug Resist*: doi:10.1089/mdr.2018.0153.
 26. Lescat M, Poirel L, Tinguely C, Nordmann P. 2019. A resazurin reductio-based assay for rapid detection of polymyxin resistance in *Acinetobacter baumannii* and *Pseudomonas aeruginosa*. *J Clin Microbiol* 57:doi:10.1128/JCM. 01563-18.
 27. Lescat M, Poirel L, Nordmann P. 2018. Rapid multiplex for detection of mcr-1 to -5 genes. *Diagn Microbiol Infect Dis* doi 10.1016/j.procbio.2018.02.003.
 28. Li J, Nation RL, Turnidge JD, Milne RW, Coulthard K, Rayner CR, Paterson DL. 2006. Colistin: the re-emerging antibiotic for multidrug-resistant Gram-negative bacterial infections. *Lancet Infect Dis* 6:589 – 601.
 29. Liu YY, Wang Y, Walsh TR, Yi LX, Zhang R, Spencer J, Doi Y, Tian G, Dong B, Huang X, Yu LF, Gu D, Ren H, Chen X, Lv L, He D, Zhou H, Liang Z, Liu JH, Shen J. 2016. Emergence of plasmid-mediated colistin resistance mechanism MCR-1 in animals and human beings in China: a microbiological and molecular biological study. *Lancet Infect Dis* 16:161–168.
 30. Lo-Ten-Foe JR, de Smet AM, Diederens BM, Kluytmans JA, van Keulen PH. 2007. Comparative evaluation of the VITEK 2, disk diffusion, Etest, broth microdilution, and agar dilution susceptibility testing methods for colistin in clinical isolates, including heteroresistant *Enterobacter cloacae* and *Acinetobacter baumannii* strains. *Antimicrob Agents Chemother* 51:3726–3730.
 31. Mali E, Florou Z, Tsilipounidaki K, Voulgaridi I, Stefos A, Xitsas S, Papagiannitsis CC, Petinaki E. 2018. Evaluation of rapid polymyxin NP to detect colistin-resistant *Klebsiella pneumoniae*. *J Microbiol Methods* 153:35-39.
 32. Matusschek E, Ahman J, Webster C, Kahlmeter G. 2017. Antimicrobial susceptibility testing of colistin- evaluation of seven commercial MCI products against standard broth microdilution for *Escherichia coli*, *Klebsiella pneumoniae*, *Pseudomonas aeruginosa* and *Acinetobacter* spp. *Clin Microbiol Infect*, doi 10.1016/j.cmi.2017.11.020.
 33. Nordmann P, Jayol A, Poirel L. 2016. Rapid detection of polymyxin resistance in *Enterobacteriaceae*. *Emerg Infect Dis* 22:1038–1043.
 34. Nordmann P, Jayol A, Poirel L. 2016. A universal culture medium for screening polymyxin-resistant Gram-negative isolates. *J Clin Microbiol* 54:1395–1399.
 35. Nordmann P, Lienhard R, Kieffer N, Clerc O, Poirel L. 2016. Plasmid-mediated colistin-resistant *Escherichia coli* in bacteremia in Switzerland. *Clin Infect Dis* 62:1322–1323.
 36. Olaitan AO, Diene SM, Kempf M, Berrazeg M, Bakour S, Gupta SK, Thongmalayvong B, Akkhavong K, Somphavong S, Paboriboune P, Chaisiri K, Komalamisra C, Adelowo OO, Fagade OE, Banjo OA, Oke AJ, Adler A, Assous MV, Morand S, Raoult D, Rolain JM. 2014. Worldwide emergence of colistin resistance in *Klebsiella pneumoniae* from healthy humans and patients in Lao PDR, Thailand, Israel, Nigeria and France owing to inactivation of the PhoP/PhoQ regulator *mgrB*: an epidemiological and molecular study. *Int J Antimicrob Agents* 44:500–507.
 37. Olaitan AO, Rolain JM. 2015. Interruption of *mgrB* in the mediation of colistin resistance in *Klebsiella oxytoca*. *Int J Antimicrob Agents* 46: 354 –355.

38. Osei Sekere J. 2019. Mcr colistin resistance gene :a systematic review of current diagnostics and detection methods. *MicrobiologyOpen* doi:10.1002/mbo3.682.
39. Pillonel T, Nordmann P, Bertelli C, Prod'Hom G, Poirel L, Greub G. Resistosome analysis of a carbapenemase (OXA-48) producing and colistin-resistant *Klebsiella pneumoniae* strain. *Antimicrob Agents Chemother* doi/10.1128/AAC.00076-18
40. Poirel L, Larpion Y, Dobias J, Stephan R, Decoussr JW, Madec JY, Nordmann P. 2017. Rapid polymyxin NP test for the detection of polymyxin resistance mediated by the *mcr-1/mcr-2* genes. *Diagn Microbiol Infect Dis*, 90:7-10.
41. Poirel L, Jayol A, Bontron S, Villegas MV, Ozdamar M, Turkoglu S, Nordmann P. 2015. The *mgrB* gene as a key target for acquired resistance to colistin in *Klebsiella pneumoniae*. *J Antimicrob Chemother* 70:75– 80.
42. Poirel L., Jayol A., Nordmann P. 2017. Polymyxins : Antimicrobial activity, susceptibility testing and resistance mechanisms encoded by plasmids and chromosome. *Clin Microb Rev* 30:557-596.
43. Poirel L, Kieffer N, Liassine N, Thanh D, Nordmann P. 2016. Plasmid- mediated carbapenem and colistin resistance in a clinical isolate of *Escherichia coli*. *Lancet Infect Dis* 16:281.
44. Rebelo AR, Bortoloia V, Kjeldgaard JS, Pedersen SK, Leekitcharoenphon P, Hansen IM, Hendriksend RS. 2018. Multiplex PCR for detection of plasmid-mediated *mcr-4* and *mcr-5* for surveillance purposes. *Eurosurveillance Weekly*, doi 10.2807/1560-7917
45. Stojanoski V, Sankaran B, Prasad BV, Poirel L, Nordmann P, Palzkill T. 2016. Structure of the catalytic domain of the colistin resistance enzyme MCR-1. *BMC Biol* 14:81.
46. Shen Z, Wang Y, Shen Y, Shen J, Wu C. 2016. Early emergence of *mcr-1* in *Escherichia coli* from food-producing animals. *Lancet Infect Dis* 16: 293.
47. Tan TY, Ng LS. 2006. Comparison of three standardized disc susceptibility testing methods for colistin. *J Antimicrob Chemother* 58:864–867.
48. Tan TY, Ng SY. 2007. Comparison of Etest, Vitek and agar dilution for susceptibility testing of colistin. *Clin Microbiol Infect* 13:541–544.
49. Xavier BB, Lammens C, Ruhel R, Kumar-Singh S, Butaye P, Goossens H, Malhotra-Kumar S. 2016. Identification of a novel plasmid-mediated colistin-resistance gene, *mcr-2*, in *Escherichia coli*, Belgium, June 2016. *Euro Surveill* 21:30280.
50. Wright MS, Suzuki Y, Jones MB, Marshall SH, Rudin SD, van Duin D, Kaye K, Jacobs MR, Bonomo RA, Adams MD. 2015. Genomic and transcriptomic analyses of colistin-resistant clinical isolates of *Klebsiella pneumoniae* reveal multiple pathways of resistance. *Antimicrob Agents Chemother* 59:536 –543.
51. Zurfuh K, Poirel L, Nordmann P, Nuesch-Inderbinen M, Hachler H, Stephan R. 2016. Occurrence of the plasmid-borne *mcr-1* colistin resistance gene in ESBL-producing *Enterobacteriaceae* in river water and imported vegetable samples in Switzerland. *Antimicrob Agents Chemother* 60:2594 –2595.